

AUTO DISABLE SYRINGE

+ Feature & Benefits

- Single handed operation and activation;
- Fingers stay behind the needle at all times;
- No change in injection technique;
- Luer slip fits into all standard hypodermic needles;
- Meets ISO standards of syringes with re-use prevention feature.

Syringe size						
1ml	2-3ml	5ml	10ml	20ml	30ml	50ml
Needle gauge						
18G	19G	20G	21G	22G	23G	
24G	25G	26G	27G	29G	30G	

Also available with other needle size on request

Unique auto-destruct function
Exclusive patent protection
WHO PQS certificates



+ How to use

<p>Draw up medicine using standard procedure</p> 	<p>Take off protector and give injection using aseptic techniques</p> 
<p>Fully depress the plunger to activate auto-destruct mechanism</p> 	<p>Dispose of syringe in sharps container</p> 